MS Tips

Mass Spectrometry Application Group Mass Spectrometry Business Unit JEOL Ltd.

No.151

JMS-T100GCV Application Data

Analysis of polymer additives using DIP-EI and FD

[Introduction]

This report shows the result of the ionization of phenolic polymer antioxidants analyzed by electron ionization (EI) with direct insertion probe(DIP) and field desorption (FD).

[Method]

Samples

Phenolic antioxidant

AO-80 (1mg/mL in methanol)

Irganox1098 (1mg/mL in methanol)

Irganox1010 (1mg/mL in methanol)

Measurement conditions

See table 1.


Fig.1 Structural formula of sample compounds

Table 1 Measurement conditions.

Instrument	JMS-T100GCV (JEOL Ltd.)
Direct probe	DIP
lonization mode	EI+ (70eV, 300μA)
Probe condition	50°C → 64°C/min → 400°C
lon source temp.	280°C
<i>m⁄z</i> range	<i>m</i> /z 50-1500
Spectrum recording time	1.0sec
Direct probe	FDP
lonization mode	FD+ (Cathode volt.: -10kV)
Probe condition	0mA → 51.2mA/min → 40mA
Ion source temp.	Heater OFF
<i>m⁄z</i> range	<i>m</i> /z 50-1500
Spectrum recording time	0.5sec



(a) DIP(<u>Direct Insertion Probe</u>)



(b) FDP(<u>Field Desorption Probe</u>)

Fig.2 Picture of direct probes.

Copyright © 2009 JEOL Ltd.




Fig.3: Mass spectra of phenolic antioxidant by DIP-EI and FD.

EI mass spectra: (a) AO-80, (b) Irganox1098, (c) Irganox1010

FD mass spectra: (A) AO-80, (B) Irganox1098, (C) Irganox1010

[Result and discussion]

The analysis time was about 5 minutes by DIP-EI and about 1 minute by FD when measurement condition in Table 1 was used.

The observed mass spectra are shown in Fig. 3. The molecular ion was observed in all mass spectra. However, the ion intensity of the molecular ion of Irganox1010 by using DIP-EI was remarkably lower than the ion intensity of other fragment ions. On the other hand, the molecular ion was observed as base peak in all mass spectra by using FD.

The JMS-T100GCV is capable to analyze antioxidants simply and quickly by using a direct injection method. This report shows that the JMS-T100GCV can be used for classic GC/MS but also for DI/MS as very useful tool for qualitative analysis.