

Thermal-decomposition Process Analysis of Sodium Formate Using TG-MS

Product used: Mass Spectrometer (MS)

Overview

Thermogravimetry (TG) is used to measure weight changes of samples under programmed heat conditions. A system combining thermogravimetry/differential thermal analysis (TG/DTA) with mass spectrometry (MS) can be used for both qualitative and quantitative analysis of gases evolved from the TG furnace into the mass spectrometer. In this application note, we show qualitative analysis of the thermal-decomposition process for sodium formate using the “STA2500 Regulus” TG system (NETZSCH) and the gas chromatography–quadrupole mass spectrometry (GC/QMS) “JMS-Q1500GC” system (JEOL).

Equipment

A 20 mg sample of sodium formate was weighed with an electronic balance and measured by TG-MS. Table 1 shows the TG-MS measurement conditions. Ionization energy was set to 17 eV, which ionized hydrogen, but not the helium used as the atmospheric gas.

Table 1. TG-MS measurement conditions.

TG conditions		MS conditions	
Furnace temp.	60°C → 20°C/min → 1000°C	Ion source temp.	250°C
Transfer line temp.	350°C	Interface temp.	300°C
Atmospheric gas flow	He, 100 mL/min	Ionization mode	EI, 17 eV
Split ratio	100:1	Ionization current	30 μA
GC conditions		Relative EM voltage	+200 V
Oven temp.	350°C	Measurement mode	Scan
Column	Blank tube 5 m × 0.25 mm i.d.	Scan range	m/z 2–100

Results

TG-MS measurement results for sodium formate are shown in Figure 1. The TG curve (= weight change), differential thermal-analysis (DTA) curve, and total ion current chromatogram (TICC) are respectively shown as green, blue, and red lines. In the DTA curve, the endothermic reaction was detected at 260°C. This was due to the melting of sodium formate (its melting point is 253°C). In the TG curve, weight loss was detected at 430°C and 570°C, and in the TICC, evolved gas peaks were detected at corresponding temperatures.




Figure 1. TG, DTA, and TICC curves of sodium formate.

Figure 2 shows the TICC and extracted-ion chromatograms (EICs) of m/z 2, 28, and 44.




Figure 2. TICC and EICs of sodium formate.

Figure 3 shows mass spectra of peaks A and B in the TICC. The mass spectrum of peak A shows m/z 2, 28, and 44, which were H_2 , CO, and CO_2 , respectively. The mass spectrum of Peak B shows m/z 28 and 44; m/z 2 was not detected.





Figure 3. Mass spectra of peaks A and B.

Consequently, TG-MS measurement results suggested the following thermal-decomposition reaction scheme [1]:



The above results suggest that TG-MS can be used to analyze chemical reactions in the thermal-decomposition process.

Reference

[1] Shunsuke Shishido and Yoshio Masuda, *Journal of the Chemical Society of Japan*, No. 1, pp. 66–70 (1976).

Copyright © 2023 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

